

Product Datasheet

BAX peptide (orb339373)

Description	BAX peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BAX
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338979 (BAX antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q07813 , Q63690 , Q07812
Entrez	581 , 12028
Expiration Date	12 months from date of receipt.