

Product Datasheet

HER2 peptide (orb339379)

Description	HER2 peptide
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ERBB2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338985 (HER2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70424, P04626, P06494
Entrez	13866, 2064
Expiration Date	12 months from date of receipt.