

Product Datasheet

CD127 peptide (orb339385)

Description	The peptide is used to CD127 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	IL7R
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338991 (CD127 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P16872, P16871
Entrez	16197, 3575
Expiration Date	12 months from date of receipt.