

Product Datasheet

Topoisomerase 2 alpha peptide (orb339390)

Description	The peptide is used to Topoisomerase 2 alpha antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	TOP2A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb338996 (Topoisomerase 2 alpha antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P11388
Entrez	7153
Expiration Date	12 months from date of receipt.