

Product Datasheet

LATS1/2 peptide (orb339397)

Description	The peptide is used to LATS1/2 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	LATS1; LATS2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339003 (LATS1/2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7TSJ6 , O95835 , Q9NRM7 , Q8BYR2
Entrez	16798 , 50523 , 9113 , 26524
Expiration Date	12 months from date of receipt.