

Product Datasheet

MEKK1 peptide (orb339398)

Description	The peptide is used to MEKK1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rabbit, Rat
Tested Applications	Blocking
Target	MAP3K1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339004 (MEKK1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q13233, Q62925, P53349
Entrez	116667, 4214
Expiration Date	12 months from date of receipt.