

Product Datasheet

Claudin 7 peptide (orb339402)

Description	The peptide is used to Claudin 7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CLDN7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339008 (Claudin 7 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9Z1L1, Q9Z261, O95471
Entrez	53624, 1366, 65132
Expiration Date	12 months from date of receipt.