

Product Datasheet

SREBP1 peptide (orb339406)

Description	The peptide is used to SREBP1 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SREBF1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339012 (SREBP1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P56720, P36956, Q9WTN3
Entrez	20787, 78968, 6720
Expiration Date	12 months from date of receipt.