

Product Datasheet

CRMP2 peptide (orb339407)

Description	The peptide is used to CRMP2 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Primate, Rat
Tested Applications	Blocking
Target	DPYSL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339013 (CRMP2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q16555, O08553, P47942
Entrez	1808, 12934, 25416
Expiration Date	12 months from date of receipt.