

Product Datasheet

IKB alpha peptide (orb339409)

Description	The peptide is used to IKB alpha antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	NFKBIA
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339015 (IKB alpha antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q63746, P25963, Q9Z1E3
Entrez	4792, 25493, 18035
Expiration Date	12 months from date of receipt.