

Product Datasheet

GRK3 peptide (orb339410)

Description	The peptide is used to GRK3 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	ADRBK2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339016 (GRK3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35626, P26819, Q3UYH7
Entrez	320129, 25372, 157
Expiration Date	12 months from date of receipt.