

Product Datasheet

CYB561D1 peptide (orb339416)

Description	The peptide is used to CYB561D1 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	CYB561D1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339022 (CYB561D1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	A2AE42, Q8N8Q1
Entrez	72023, 284613
Expiration Date	12 months from date of receipt.