

Product Datasheet

ACTN1 peptide (orb339417)

Description	The peptide is used to ACTN1 antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	ACTN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339023 (ACTN1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7TPR4, P12814, Q9Z1P2
Entrez	87, 109711, 81634
Expiration Date	12 months from date of receipt.