

Product Datasheet

RPS20 peptide (orb339423)

Description	The peptide is used to RPS20 antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	RPS20
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339029 (RPS20 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P60867, P60866, P60868
Entrez	100362149, 6224, 122772, 67427, 100359951, 100362684
Expiration Date	12 months from date of receipt.