

Product Datasheet

ATP5F1 peptide (orb339424)

Description	The peptide is used to ATP5F1 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ATP5F1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339030 (ATP5F1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P24539, Q9CQQ7
Entrez	515, 11950
Expiration Date	12 months from date of receipt.