

Product Datasheet

B3GALTL peptide (orb339425)

Description	The peptide is used to B3GALTL antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	B3GALTL
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339031 (B3GALTL antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q6Y288
Entrez	145173
Expiration Date	12 months from date of receipt.