

## Product Datasheet

### CSGALNACT2 peptide (orb339429)

|                            |                                                                                                                                                                        |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to CSGALNACT2 antibody reactivity.                                                                                                                 |
| <b>Reactivity</b>          | Human                                                                                                                                                                  |
| <b>Tested Applications</b> | Blocking                                                                                                                                                               |
| <b>Target</b>              | CSGALNACT2                                                                                                                                                             |
| <b>Form/Appearance</b>     | Lyophilized powder                                                                                                                                                     |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                                                                                                                           |
| <b>Note</b>                | For research use only                                                                                                                                                  |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb339035 (CSGALNACT2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. |
| <b>Purity</b>              | > 85%                                                                                                                                                                  |
| <b>Source</b>              | Synthetic                                                                                                                                                              |
| <b>Uniprot ID</b>          | <b>Q8N6G5</b>                                                                                                                                                          |
| <b>Entrez</b>              | <b>55454</b>                                                                                                                                                           |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                        |