

Product Datasheet

NEK7 peptide (orb339444)

Description	The peptide is used to NEK7 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	NEK7
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339050 (NEK7 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8TDX7, D3ZBE5, Q9ES74
Entrez	360850, 140609, 59125
Expiration Date	12 months from date of receipt.