

Product Datasheet

HNF1 beta peptide (orb339450)

Description	The peptide is used to HNF1 beta antibody reactivity.
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	HNF1B
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339056 (HNF1 beta antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35680, P23899, P27889
Entrez	6928, 25640, 21410
Expiration Date	12 months from date of receipt.