

## Product Datasheet

### Cortactin (AcK235) peptide (orb339452)

|                            |                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The peptide is used to Cortactin (AcK235) antibody reactivity.                                                                                                                 |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                                                              |
| <b>Tested Applications</b> | Blocking                                                                                                                                                                       |
| <b>Target</b>              | CTTN                                                                                                                                                                           |
| <b>Form/Appearance</b>     | Lyophilized powder                                                                                                                                                             |
| <b>Storage</b>             | Shipped at 4°C. Store at -20°C for one year.                                                                                                                                   |
| <b>Note</b>                | For research use only                                                                                                                                                          |
| <b>Application notes</b>   | This peptide was used as immunogen to produce orb339058 (Cortactin (AcK235) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. |
| <b>Purity</b>              | > 85%                                                                                                                                                                          |
| <b>Source</b>              | Synthetic                                                                                                                                                                      |
| <b>Uniprot ID</b>          | <b>Q60598, Q14247</b>                                                                                                                                                          |
| <b>Entrez</b>              | <b>13043, 2017</b>                                                                                                                                                             |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                |