

Product Datasheet

Cortactin peptide (orb339453)

Description	The peptide is used to Cortactin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	CTTN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339059 (Cortactin antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q60598, Q14247
Entrez	13043, 2017
Expiration Date	12 months from date of receipt.