

Product Datasheet

XPO1 (AcK568) peptide (orb339454)

Description	The peptide is used to XPO1 (AcK568) antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	XPO1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339060 (XPO1 (AcK568) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O14980, Q80U96, Q6P5F9
Entrez	7514, 85252, 103573
Expiration Date	12 months from date of receipt.