

Product Datasheet

14-3-3-pan (AcK51/49) peptide (orb339455)

Description	The peptide is used to 14-3-3-pan (AcK51/49) antibody reactivity.
Reactivity	Human, Mouse, Rat, Sheep, Zebrafish
Tested Applications	Blocking
Target	YWHAB; YWHAG; YWHAQ; YWHAZ; SFN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339061 (14-3-3-pan (AcK51/49) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61982, P31946, P68254, P68255, P27348, P35213, P63101, P63102, P31947, P61981, P63104, Q9CQV8, P61983
Entrez	22630, 56011, 56010, 7534, 22631, 54401, 25578, 7529, 25577, 22628, 7532, 2810, 10971
Expiration Date	12 months from date of receipt.