

Product Datasheet

p53 peptide (orb339469)

Description	The peptide is used to p53 antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	TP53
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339075 (p53 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P10361, P04637
Entrez	7157, 24842
Expiration Date	12 months from date of receipt.