

Product Datasheet

Alpha-tubulin (AcK352) peptide (orb339477)

Description	The peptide is used to Alpha-tubulin (AcK352) antibody reactivity.
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TUBA1A; TUBA1B; TUBA1C; TUBA3C; TUBA4A; TUBA8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339083 (Alpha-tubulin (AcK352) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry. Directions for Use: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P05213, Q6AYZ1, Q13748, P68370, Q6P9V9, Q9BQE3, Q9NY65, P68363, Q71U36, P68369, P68366, P68373
Entrez	1962, 64158, 51807, 7278, 113457, 10376, 7277, 500929, 22146, 22142, 84790, 22143, 7846, 100909441
Expiration Date	12 months from date of receipt.