

Product Datasheet

FGF12 peptide (orb339487)

Description	The peptide is used to FGF12 antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FGF12
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339093 (FGF12 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P61328, P61150, P61329
Entrez	2257, 14167, 170630
Expiration Date	12 months from date of receipt.