

Product Datasheet

Apelin peptide (orb339492)

Description	The peptide is used to Apelin antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	APLN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339098 (Apelin antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9R0R3, Q9ULZ1, Q9R0R4
Entrez	58812, 8862, 30878
Expiration Date	12 months from date of receipt.