

Product Datasheet

Angiopoietin 2 peptide (orb339497)

Description	Angiopoietin 2 peptide
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	ANGPTL1; ANGPTL2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339103 (ANGPTL1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q95841, Q9UKU9, Q640P2, Q9R045
Entrez	72713, 26360, 9068, 23452
Expiration Date	12 months from date of receipt.