

Product Datasheet

STRA8 peptide (orb339501)

Description	The peptide is used to STRA8 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	STRA8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339107 (STRA8 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q7Z7C7, P70278
Entrez	20899, 346673
Expiration Date	12 months from date of receipt.