

Product Datasheet

CCL11 peptide (orb339509)

Description	The peptide is used to CCL11 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	CCL11
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339115 (CCL11 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P51671, P48298
Entrez	6356, 20292
Expiration Date	12 months from date of receipt.