

Product Datasheet

TGF beta 4 peptide (orb339543)

Description	The peptide is used to TGF beta 4 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	LEFTY1; LEFTY2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339149 (TGF beta 4 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	O75610 , O00292
Entrez	7044 , 10637
Expiration Date	12 months from date of receipt.