

Product Datasheet

Neurturin peptide (orb339545)

Description	The peptide is used to Neurturin antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	NRTN
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339151 (Neurturin antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P97463, Q99748
Entrez	4902, 18188
Expiration Date	12 months from date of receipt.