

Product Datasheet

BMP3A peptide (orb339546)

Description	The peptide is used to BMP3A antibody reactivity.
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	BMP3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339152 (BMP3A antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8BHE5, P12645, P49002
Entrez	25667, 110075, 651
Expiration Date	12 months from date of receipt.