

Product Datasheet

CCL4L1 peptide (orb339549)

Description	The peptide is used to CCL4L1 antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	CCL4L1; CCL4L2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339155 (CCL4L1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8NHW4
Entrez	9560, 388372
Expiration Date	12 months from date of receipt.