

Product Datasheet

S100A1 peptide (orb339558)

Description	S100A1 peptide
Reactivity	Human
Tested Applications	Blocking
Target	S100A1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339164 (S100-A1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P23297
Entrez	6271
Expiration Date	12 months from date of receipt.