

Product Datasheet

TP53INP2 peptide (orb339559)

Description	The peptide is used to TP53INP2 antibody reactivity.
Reactivity	Human, Mouse
Tested Applications	Blocking
Target	TP53INP2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339165 (TP53INP2 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q8IXH6, Q8CFU8
Entrez	58476, 68728
Expiration Date	12 months from date of receipt.