

Product Datasheet

Histamine H1 Receptor peptide (orb339567)

Description	The peptide is used to Histamine H1 Receptor antibody reactivity.
Reactivity	Human
Tested Applications	Blocking
Target	HRH1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339173 (Histamine H1 Receptor antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P35367
Entrez	3269
Expiration Date	12 months from date of receipt.