

Product Datasheet

MC5 Receptor peptide (orb339568)

Description	The peptide is used to MC5 Receptor antibody reactivity.
Reactivity	Human, Mouse, Sheep, Zebrafish
Tested Applications	Blocking
Target	MC5R
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339174 (MC5 Receptor antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P33032, P41149
Entrez	4161
Expiration Date	12 months from date of receipt.