

Product Datasheet

Neurokinin B receptor Peptides (orb339569)

Description	Neurokinin B receptor Peptides
Reactivity	Human
Tested Applications	Blocking
Target	TACR3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339175 (Tachykinin Receptor 3 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P29371
Entrez	6870
Expiration Date	12 months from date of receipt.