

Product Datasheet

VAV1 peptide (orb339582)

Description	The peptide is used to VAV1 antibody reactivity.
Reactivity	Human, Rat
Tested Applications	Blocking
Target	VAV1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339188 (VAV1 antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	P15498, P54100
Entrez	7409
Expiration Date	12 months from date of receipt.