

Product Datasheet

eNOS (phospho-T494) peptide (orb339587)

Description	The peptide is used to eNOS (pT494) antibody reactivity.
Reactivity	Human, Mouse, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	NOS3
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb339193 (eNOS (phospho-T494) antibody)The quality of the peptide was evaluated by reversed-phase HPLC and by mass spectrometry.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q62600, P29474, P70313
Entrez	24600, 4846, 18127
Expiration Date	12 months from date of receipt.