



Product Datasheet

ORM1 antibody (orb339607)


Description	ORM1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human ORM1 protein. The protein is available as (orb315543) was used as immunogen.
Target	ORM1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	24 kDa
Dilution Range	IHC-P: 1:100-1000, WB: 1:100-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



WB analysis of rat small intestines (lanes 1-6)



IHC-P staining of mouse skin tissue using ORM1 antibody