



Product Datasheet

alpha 1 Antitrypsin antibody (orb339609)


Description	alpha 1 Antitrypsin antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant human alpha 1 Antitrypsin protein. The protein is available as (orb315533) was used as immunogen.
Target	alpha 1 Antitrypsin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	42.8 kDa
Dilution Range	IHC-P: 1:100-1000, WB: 1:100-2000
Expiration Date	12 months from date of receipt.



IHC-P
staining of
mouse skin
tissue usin...



WB analysis
of human
lung cancer
(lane 1...



IHC-P
staining of
rat ovary
tissue
using...