



## Product Datasheet

### PEDF antibody (orb339612)

|                            |                                                                                                                |
|----------------------------|----------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | PEDF antibody                                                                                                  |
| <b>Species/Host</b>        | Rabbit                                                                                                         |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                   |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                      |
| <b>Immunogen</b>           | Recombinant human PEDF protein. The protein is available as (orb315544) was used as immunogen.                 |
| <b>Target</b>              | PEDF                                                                                                           |
| <b>Form/Appearance</b>     | 10 mM PBS, 0.02% sodium azide                                                                                  |
| <b>Concentration</b>       | - 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml                                                 |
| <b>Storage</b>             | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| <b>Note</b>                | For research use only                                                                                          |
| <b>Isotype</b>             | IgG                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                     |
| <b>Purity</b>              | Polyclonal antibodies are purified by peptide affinity chromatography                                          |
| <b>MW</b>                  | 46 kDa                                                                                                         |
| <b>Dilution Range</b>      | WB: 1:100:2000, IHC-P: 1:50-500                                                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                |



Western blot analysis of PEDF protein 50...



Immunohistochemical staining of rat ovar...



WB analysis of human lung cancer (lane 1...