

## Product Datasheet

### Human dsDNA ELISA Kit (orb339621)

|                          |                                                                 |
|--------------------------|-----------------------------------------------------------------|
| <b>Description</b>       | Human anti-double stranded DNA (dsDNA) antibody (IgG) ELISA Kit |
| <b>Reactivity</b>        | Human                                                           |
| <b>Range</b>             | Request Information                                             |
| <b>Target</b>            | dsDNA Ab (IgG)                                                  |
| <b>Note</b>              | For research use only                                           |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                    |
| <b>Sample Types</b>      | serum                                                           |
| <b>Sensitivity</b>       | Request Information                                             |
| <b>Expiration Date</b>   | Please enquire.                                                 |