

## Product Datasheet

### Human HCG ELISA Kit (orb339650)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Description</b>       | Human HCG ELISA Kit                           |
| <b>Reactivity</b>        | Human                                         |
| <b>Range</b>             | 8 mIU/mL-240 mIU/mL                           |
| <b>Target</b>            | hCG                                           |
| <b>Note</b>              | For research use only                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                  |
| <b>Sample Types</b>      | serum, plasma, urine, cell culture supernates |
| <b>Uniprot ID</b>        | <b>PODN86</b>                                 |
| <b>Sensitivity</b>       | 2 mIU/mL                                      |
| <b>Expiration Date</b>   | Please enquire.                               |