

Product Datasheet

Human MMP-9 ELISA Kit (orb339657)

Description	Human matrix metalloproteinase 9/Gelatinase B, MMP-9 ELISA Kit
Reactivity	Human
Range	0.312 ng/mL-20 ng/mL
Target	MMP9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates, urine
Uniprot ID	P14780
Sensitivity	0.284 ng/mL
Expiration Date	Please enquire.