

Product Datasheet

Mouse apo-A5 ELISA Kit (orb339666)

Description	Mouse Apolipoprotein A5 (apo-A5) ELISA Kit
Reactivity	Mouse
Range	31.2 pg/mL-2000 pg/mL
Target	APOA5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8C7G5
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.