

Product Datasheet

Mouse MMP-9 ELISA Kit (orb339667)

Description	Mouse matrix metalloproteinase 9/Gelatinase B, MMP-9 ELISA Kit
Reactivity	Mouse
Range	0.625 ng/mL-40 ng/mL
Target	MMP9
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P41245
Sensitivity	0.552 ng/mL
Expiration Date	Please enquire.