

Product Datasheet

Mouse VEGF ELISA Kit (orb339669)

Description	Mouse Vascular Endothelial cell Growth Factor, VEGF ELISA KIT
Reactivity	Mouse
Range	3.906 pg/mL-250 pg/mL
Target	VEGFA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, cell culture supernates, tissue homogenates
Uniprot ID	Q00731
Sensitivity	0.857 pg/mL
Expiration Date	Please enquire.